BRAUN

ExactFit™ 1

BUA5000

Braun ExactFit™ 1 BUA5000

Upper arm blood pressure monitor

GB	English.....	2
DE	Deutsch	11
DK	Dansk.....	21
ES	Español	30
FI	Suomi.....	39
FR	Français.....	48
IT	Italiano.....	57
NL	Nederlands.....	66
NO	Norsk.....	75
PL	Polski.....	84
SA	اللغة العربية.....	101
SE	Svenska.....	102
SI	Slovenščina	111
TR	Türkçe.....	120

Intended use of Braun ExactFit™ 1

Braun ExactFit™ 1 is a blood pressure monitor that has been developed for easy and accurate blood pressure measurements from the upper arm. The measuring accuracy of Braun ExactFit™ 1 was tested at the time of manufacture and was proven by clinical research in accordance with ESH.¹

This device has not been tested for and hence is not intended for neonatal and infant subjects, pregnant, including pre-eclamptic users.

What you should know about blood pressure

Blood pressure constantly changes throughout the day. It rises sharply in the early morning and declines during the late morning. Blood pressure rises again in the afternoon and finally drops to a low level at night. Also, it may vary in a short period of time. Therefore, readings from successive measurements can fluctuate.

This device will display your blood pressure reading referring to two values: systolic and diastolic. The systolic blood pressure (the upper number) indicates how much pressure your blood is exerting against your artery walls when the heart beats. The diastolic blood pressure (the lower number) indicates how much pressure your blood is exerting against your artery walls while the heart is resting between beats. The pulse value is also displayed after each measurement.

Blood pressure measured in a doctor's office only provides a momentary value. Repeated measurements at home better reflect one's actual blood pressure values under everyday conditions.

Moreover, many people have a different blood pressure when they measure at home, because they tend to be more relaxed than when in the doctor's office. Regular blood pressure measurements taken at home can provide your doctor with valuable information on your normal blood pressure values under actual «everyday» conditions.

The ESH (European Society of Hypertension) and WHO (World Health Organisation) have set up the following standard blood pressure values when measured with a resting pulse at home:

Blood pressure (mmHg)	Normal values	Grade 1 mild hypertension	Grade 2 moderate hypertension	Grade 3 severe hypertension
SYS = systole (upper value)	up to 134	135-159	160-179	≥180
DIA = diastole (lower value)	up to 84	85-99	100-109	≥110

¹ European Society of Hypertension

Safety information and important precautions

To ensure accurate measuring results, read the complete instructions for use carefully.

People suffering from cardiac arrhythmia, vascular constriction, arteriosclerosis in the extremities, diabetes or users of cardiac pacemakers should consult their doctor before measuring their blood pressure themselves. Deviations in blood pressure values may occur in such cases.

If your blood pressure reading is in the severe hypertension range with values of 180 mmHg systolic and above, and/or 110 mmHg diastolic and above; and you find you have one or two readings in this range - please see your doctor straight away.

If you are under medical treatment or taking any medication, please consult your doctor first.

The use of this blood pressure monitor is not intended as a substitute for a consultation with your doctor.

Avoid compression or restriction of the connection tubing. Continuous cuff pressure due to the connection tubing kinking can affect blood flow and cause possible injury to the user.

It is not recommended to use this device on the arm of the same side of the body on which the user had a mastectomy.

Measuring blood pressure too frequently can cause injury to the user due to blood flow interference.

It is not recommended to apply the cuff to a limb on which there is intravascular access or on which an arterio-venous (A-V) shunt is present. This can cause temporary interference of blood flow and potential injury to the user.

DO NOT apply the cuff over a wound as this can cause further injury.


It is not recommended to use the device while simultaneously using other medical electrical (ME) equipment on the same limb.

Product description

- 1 LCD display
- 2 Power (start/stop) button 
- 3 Hose port
- 4 Battery compartment
- 5 Connector
- 6 Arm cuff
- 7 Air hose

Inserting batteries

- Use alkaline batteries, type AA 1.5V as supplied with this product.
- Remove the battery compartment cover and insert four batteries with correct polarity (see symbol in the battery compartment).

-  Batteries should not be disposed of in the household waste, but at appropriate collection sites or at your retailer.

Key rules for accurate blood pressure measurement

Always take readings at the same time of day, under the same conditions.

Do not take a measurement within 30 minutes of smoking, drinking coffee or tea or any form of exertion. These factors will influence the measurement.

Ensure at least 5 minutes of quiet rest before measurements. If possible, it is recommended to rest 10 minutes before taking the measurement. Try to always use the same resting time before taking a measurement.

Always measure on the same arm. The left arm is recommended. Do not move or speak. Sit comfortably, keep your legs uncrossed and feet flat on the floor, back and arm supported.

Place the cuff in direct contact with the skin. Do not wrap the cuff over a jacket, sweater or shirt.

If you want to take multiple measurements, wait approximately 3 minutes before repeating a measurement.

Applying the arm cuff


Fig. 1


Fig. 2


Fig. 3

- Plug the air hose into the connector (Fig. 1).
- The cuff is located correctly when the velcro is on the outside of the cuff and the metal ring does not touch the skin (Fig. 2).

- Put your arm through the cuff loop. The bottom of the cuff should be approximately 2-3 cm above elbow (Fig. 3). The tube should lie over the brachial artery on the inside of the arm (Fig. 4).


Fig. 4


Fig. 5


Fig. 6

- Pull the cuff so that the top and bottom edges are tightened around your arm.
- When the cuff is positioned properly, press the Velcro firmly against the pile side of the cuff.
- This cuff is suitable for use if the arrow heads fall within the stripe pattern on the other side of the cuff when the cuff is tightened around your arm (Fig. 5).
- Sit on a chair and place your arm on the table so that the cuff is at the same level as your heart (Fig. 6).

Taking a measurement

- Sit on a chair with your feet flat on the floor.
- Wrap the cuff snugly around your arm (see "Applying the arm cuff" section above, Fig. 3-6).
- Do not move the unit during measurement, or the proper measurement will not be achieved.
- Press power button  to begin. The measurement starts automatically. While measurement is in progress, do not move or speak.

Note: If necessary to interrupt a measurement, press the power button  at any time. The device immediately cancels the measurement, lowers cuff pressure and automatically turns off.

- The cuff will first inflate and then deflate. At the end of the measurement there are 2 final beep sounds and the reading is displayed:

Upper (systolic) value
Lower (diastolic) value
Pulse

- Press the power button  to turn the product off. Otherwise the

product will turn off automatically after 1 minute.

14. Test results will be saved into the memory automatically (see memory function below).

This device is not intended as a substitute for regular check-ups by your doctor, please continue to visit your doctor on a regular basis for a professional reading.

Memory function

Your ExactFit™ 1 blood pressure monitor will automatically store the latest measurement. Storing is done automatically after the completed measurement.

To recall the reading, press and hold the power button  for at least 4 seconds while the unit is switched off. When all the display elements are shown, release the power button and the screen will show the last reading with an **M** symbol.

Battery charge indicator


Batteries almost discharged

When the batteries are approximately 75% used, a battery warning appears. This icon blinks 3 times first when the device is turned on, then goes to measuring mode. The monitor will continue to measure reliably, but new batteries should be obtained.

Batteries discharged—replacements required

 When the batteries are fully discharged, the battery symbol will flash empty along with 3 short beep sounds as soon as the device is switched on. You cannot take any further measurement and must replace the batteries. Refer to the section called "Inserting batteries".

Storage and cleaning

- Do not expose device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight inflatable bladder. Handle this carefully and avoid all types of straining through twisting or buckling.
- Clean the device with a soft dry cloth. Stains on the cuff can be carefully removed with a damp cloth. Do not use gasoline, thinners or similar solvents. Do not wash or dry clean the cuff. Do not store the

product in a place exposed to direct sunlight or high humidity.

- Do not drop the device or treat it roughly in any way. Avoid strong vibrations.
- NEVER open the device! This will make the manufacturer's warranty invalid and could adversely affect the device's accuracy!

Calibration

This product has been calibrated at the time of manufacture. If used according to the instructions for use, periodic recalibration is not required. If at any time you question the measuring accuracy, please contact our service representative (see last page for contact information).

What to do if ...

Error	Description	Potential cause and remedy
	Batteries are low	Get ready to replace batteries soon with new ones.
	Batteries are fully discharged	Insert new batteries or check that batteries are correctly inserted.
«Err 1»	Signal too weak	The pulse signals on the cuff are too weak. Re-position the cuff and repeat the measurement.
«Err 2»	Error signal	During the measurement, error signals were detected by the cuff, caused for instance by movement or muscle tension. Repeat the measurement, keeping the arm still.
«Err 3»	No pressure in the cuff	An adequate pressure cannot be generated in the cuff. A leak may have occurred. Check that the cuff is correctly connected and is not too loose. Replace the batteries if necessary. Repeat the measurement.
«Err 5»	Abnormal result	The measuring signals are inaccurate and no result can therefore be displayed. Read through the checklist for performing reliable measurements and then repeat the measurement.
«Hl»	Pulse or cuff pressure too high	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.
«Lo»	Pulse too low	The pulse is too low (less than 40 beats per minute). Repeat the measurement.

Error	Description	Potential cause and remedy
Readings from repeated measurements differ considerably.	Blood pressure is a fluctuating value. For healthy adults, deviations of 10 to 20 mmHg are possible.	Reminder: always use the same arm for measurements!
	Cuff is not properly fitted.	Ensure that the cuff is fitted according to the "Applying the arm cuff" section.
	Readings were not consistently taken at heart level.	For each measurement, make sure that the cuff is at heart level.
	Talking, coughing, laughing, moving etc. when measuring will influence the reading.	While measuring, relax, keep still, do not move or speak.
Blood pressure values measured at your doctor's office differ from your measurements.	Doctor visits can frequently cause anxiety.	Record the daily development of the measured values and consult your doctor.
Display remains blank, or is unusual, when instrument is switched on.	Batteries are not installed properly.	Check batteries for correct polarity.
Device frequently fails to measure blood pressure values, or values measured are too low (or too high).	Cuff may not be properly positioned.	Ensure that the cuff is fitted according to the "Applying the arm cuff" section.

Specifications

Model:	BUA5000
Operation principle:	Oscillometric method
Display:	Liquid crystal display
Range of measurement:	Pressure 30~280 mmHg Pulse 40~199 beats/minute
Laboratory accuracy:	±3 mmHg (cuff pressure) ±5% of reading (pulse rate)
Clinical accuracy:	<5 mmHg systematic offset <8 mmHg standard deviation
Battery life:	400 measurements
Auto power off:	Whenever not used for 1 minute
Accessories:	4 batteries, 1 arm cuff with tube, instruction manual
Inflation:	Automatic
Power supply:	4 batteries, type AA 1.5V
Operating temperature/ atmospheric pressure/ humidity:	10 °C to 40 °C (50 °F to 104 °F) 860 to 1060 hPa (0.849 to 1.046 atm) 15% to 90% relative humidity, non-condensing
Storage temperature/ atmospheric pressure/ humidity:	-20 °C to 55 °C (-4 °F to 131 °F) 860 to 1060 hPa (0.849 to 1.046 atm) 15% to 90% relative humidity, non-condensing
Cuff:	Universal cuff (22-42 cm)
Service life:	5 years
Warranty:	2 years

If device is not used within specified temperature and humidity ranges the technical accuracy cannot be ensured.

Type BF applied parts

Operating temperature

See instructions for use

Storage temperature

Keep dry

Storage humidity/
operating humidity

Subject to change without notice.

Internally powered equipment.

Continuous operation.

This product conforms to the provisions of the EC directive 93/42/EEC (Medical Device Directive).

This product carries the CE mark and is manufactured in conformity with the RoHS Directive 2011/65/EU.

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC. For detailed description of EMC requirements visit www.helenofroy.com/emea-en/support/. Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.

Please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at your local retailer or at appropriate collection points provided in your country.

Warranty

Please read all instructions before attempting to use this device. Please retain the receipt as proof of and date of purchase. The receipt must be presented when making any claim within the relevant warranty period. Any claim under the warranty will not be valid without a proof of purchase.

Your device is warrantied for two years (2 years) from date of purchase.

This warranty covers defects in materials or workmanship that occur under normal use; defective devices meeting these criteria will be replaced free of charge.

The warranty DOES NOT cover defects or damage resulting from abuse or failure to follow the user instructions. The warranty becomes void if the device is opened, tampered with, or used with non-Braun branded parts or accessories, or if repairs are undertaken by unauthorised persons.

Accessories and consumables are excluded from any warranty.

For support requests, please visit www.helenofroy.com/emea-en/support/ or find service contact information at the end of this owner's manual.

This warranty is applicable to Europe, Middle East and Africa only.

UK Only: This does not affect your consumer statutory rights.

The LOT and SN of your device are printed on the rating label on the back of the product. The first 3 numeric digits of the LOT No. represents the day of the year of manufacture. The next 2 numeric digits represent the last two numbers of the calendar year of manufacture and the letters at the end designate the manufacturer of the product. E.g. LOT 15614onb, this product was made on day 156, year 2014, at manufacturer identifier onb.

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☎ +44 (0) 207 949 0115

www.BraunHealthcare.com

Certain trademarks used under license from
The Procter & Gamble Company or its affiliates.

ExactFit™ is a trademark of Helen of Troy
Limited and/or its affiliates.

Kaz Europe Sàrl
A Helen of Troy Company
Place Chauderon 18
CH-1003 Lausanne
Switzerland

Made and
printed in China

© 2021 All rights reserved.
P/N A005565R1
30MAR21

Brand: Braun
Category: BPM
Model: **BUA5000EUV1
Series**

Artwork Part #: A005565R1
Die Line Part #: N/A
Subject: Owner's Manual
Region: EMEA
Flat Size: MM: W180x H160

Folded Size: MM: W90 x H160

Scale: 1/1
Material: 50# gloss or matte text (inside pages)
70# gloss or matte text (cover)

Page Count: 132
Revision: 4

Date: **22APR21**
Release Date: 28APR21
Re-release Date:

Colors:

Dielines (Do not print)

Spot Colors

Black
100%PMS
Goe 63-1-4U**Special Instructions:**

Quality Requirement of Artwork and Quality
Clarification Process of Artwork Printing:
Meet Eng-QS-06&02

Helen of Troy
Creative Services
Marlborough, MA 01752, USA
+1 508 490 7000